



COSMOSIL

COSMOSIL UFLC Columns

COSMOSIL UFLC Columns with 2.5 µm particles are a powerful tool for very fast and efficient separation. Furthermore they can be used under 1/2 pressure of competitors' 2 µm columns.

Packing Material	2.5C ₁₈ -MS-II	2.5Cholester	2.5πNAP
Silica Gel	High Purity Porous Spherical Silica		
Average Particle Size	2.5 µm		
Average Pore Size	approx. 130 Å		
Specific Surface Area	approx. 330m ² /g		
Stationary Phase	 Octadecyl Group	 Cholesteryl Group	 Naphtylethyl Group
Bonding Type	Monomeric		
Main Interaction	Hydrophobic Interaction	Hydrophobic Interaction Molecular Shape Selectivity	Hydrophobic Interaction π-π Interaction
End Capping Treatment	Near-perfect Treatment		
Features	- Multi-purpose C₁₈ column - Suitable for basic compounds.	- Usable under the same condition as C₁₈. - High molecular shape selectivity	- Stronger π-π interaction than Phenyl columns.

COSMOSIL 2.5C₁₈-MS-II

Application Data


COSMOSIL 2.5C₁₈-MS-II (3.0 mm I.D. x 75 mm) enables the separation only in 2 minutes.

• Prednisone and Salicylic Acid

COSMOSIL Application Data

Column: 2.5C₁₈-MS-II
 Column size: 3.0mm I.D.-75mm
 Mobile phase: Acetonitrile/ H₂O = 30/70
 Flow rate: 1.0 ml/min
 Temperature: 40°C
 Detection: UV242nm


Sample: Prednisone (0.02mg/ml)
 Inj. Vol. 5 µl



COSMOSIL Application Data

Column: 2.5C₁₈-MS-II
 Column size: 3.0mm I.D.-75mm
 Mobile phase: Acetonitrile/ 20mmol/l Phosphate Buffer(pH2.5) = 30/70
 Flow rate: 1.0 ml/min
 Temperature: 40°C
 Detection: UV296nm


Sample: Salicylic Acid (0.33mg/ml)
 Inj. Vol. 0.5 µl




Pressure comparison with competitors' 2 μ m columns

COSMOSIL 2.5C₁₈-MS-II can be used under 1/2 pressure of competitors' 2 μ m columns.


• 3.0 mm I.D. x 75 mm



• 2.0 mm I.D. x 50 mm



The relations between flow rate, retention time, pressure and theoretical plate number



- Low back pressure (2.5 μm silica gel)
- Excellent molecular shape selectivity
- Usable under the same condition as ODS

Improved Separation


COSMOSIL 2.5Cholester offers improved resolution for compounds difficult to analyze with C_{18} without changing analytical condition.

• Catechins


COSMOSIL Application Data

Column size: 3.0mm I.D.-75mm
 Mobile phase: A: Acetonitrile/ 20mmol/l Phosphate Buffer(pH2.5) = 10/90
 B: Acetonitrile/ 20mmol/l Phosphate Buffer(pH2.5) = 30/70
 B conc. 0 $\rightarrow$ 100% 5min Linear gradient
 Mixer: 0.5ml
 Flow rate: 1.0 ml/min
 Temperature: 40°C
 Detection: UV280nm
 Sample:
 1: Gallocatechin(GC) (0.40mg/ml)
 2: Caffeine (0.04mg/ml)
 3: Epigallocatechin(EGC) (0.40mg/ml)
 4: Catechin (C) (0.20mg/ml)
 5: Epicatechin(EC) (0.20mg/ml)
 6: Epigallocatechin gallate(EGCG) (0.10mg/ml)
 7: Gallocatechin gallate(GCG) (0.20mg/ml)
 8: Epicatechin gallate(ECG) (0.10mg/ml)
 9: Catechin gallate(CG) (0.10mg/ml)
 Injection Vol. 1.0 μl

Company A 2 μm C_{18}



COSMOSIL 2.5Cholester




• Saikosaponins


COSMOSIL Application Data

Mobile phase: Acetonitrile/ 0.05% NaH_2PO_4 aq. = 3/5
 Flow rate: 0.7 ml/min
 Temperature: 50°C
 Detection: UV206nm
 Sample:
 1: Saikosaponin c
 2: Saikosaponin h
 3: Saikosaponin a
 4: Saikosaponin b2
 5: Saikosaponin b1
 6: Saikosaponin d
 Injection Vol. 1.0 μl

Company C 1.7 μm C_{18}
 (2.1mm I.D.-50mm)



COSMOSIL 2.5Cholester
 (2.0mm I.D.-50mm)



• Stilbene (Geometrical Isomers)

COSMOSIL Application Data

Column:
 Column size: 3.0mm I.D.-75mm
 Mobile phase: Acetonitrile/ H_2O = 70/30
 Flow rate: 1.0 ml/min
 Temperature: 40°C
 Detection: UV254nm
 Sample:
 1: *cis*-Stilbene (0.10mg/ml)
 2: *trans*-Stilbene (0.20mg/ml)
 Injection Vol. 1.0 μl

Company A 2 μm C_{18}



COSMOSIL 2.5Cholester




- Low back pressure (2.5 μ m silica gel)
- Naphthalene bonded stationary phase
- Stronger π - π interactions than Phenyl columns


Improved Separation

COSMOSIL 2.5 π NAP provides greater performance in separating positional isomers and other closely related compounds which are difficult to analyze with C₁₈.

• Tolunitriles (Positional Isomers)



• Tocopherols (Positional Isomers)



Ordering Information

Size	COSMOSIL 2.5C ₁₈ -MS-II	COSMOSIL 2.5Cholester	COSMOSIL 2.5 π NAP
I.D.xLength (mm)	Product Number	Product Number	Product Number
2.0x50	08994-31	09000-01	06062-91
2.0x75	08995-21	09047-11	06051-31
2.0x100	08996-11	09048-01	06052-21
3.0x50	08997-01	09049-91	06054-01
3.0x75	08998-91	09050-51	06055-91
3.0x100	08999-81	09051-41	06057-71

For research use only, not intended for diagnostic or drug use.